

HYDRON SOFTINTS™

-  SAPPHIRE
-  CHESTNUT
-  EMERALD
-  AQUA
-  AMETHYST
-  CORAL
-  AMBER

HYDRON SOFTINTS™ NOW EXCLUSIVELY FROM CLPL

- Zero 6™ Sphere Tint
- OmniFlex® Sphere Tint
- Mini™ Tint
- H67™ Tint
- Zero 6T™ Toric Tint
- OmniFlex® Toric Tint










CONTACT LENSES	PRICE PER LENS (£)	3 MONTHS (£)	6 MONTHS (£)	12 MONTHS (£)	WATER CONTENT
SPHERICAL					
Hydron OmniFlex® SOFTINTS™					70.0%
Hydron Zero 6™ SOFTINTS™					38.6%
Hydron Zero 6™ SOFTINTS™ (Rx)					38.6%
Hydron Mini™ SOFTINTS™					38.6%
Hydron H67™ SOFTINTS™					67.0%
TORIC					
Hydron OmniFlex® Toric SOFTINTS™					70.0%
Hydron Zero 6T™ Toric SOFTINTS™					38.6%

N O W e x c l u s i v e l y

EFFECTIVE FROM 3RD FEBRUARY 2003

● SAPPHIRE
 ● CHESTNUT
 ● EMERALD
 ● AQUA
 ● AMETHYST
 ● CORAL
 ● AMBER

ADDITIONAL OPTION: ● COSMETIC BLACK PUPIL

BASE CURVE (mm)	DIAMETER (mm)	SPHERICAL POWERS (D)	CYLINDRICAL POWERS (D)	AXES	SOFTINTS™
8.4, 8.8	14.3	±10.00 (0.25 steps) ±10.50 to ±20.00 (0.50 steps)	-	-	 (10%, 20%)
8.4 to 9.3 (0.3 steps)	14.0	±10.00 (0.25 steps) ±10.50 to ±20.00 (0.50 steps)	-	-	 (5%, 10%, 20%, 40%)
8.1 to 9.6 (0.3 steps)	13.5, 14.0, 14.5, 15.0	±10.00 (0.25 steps) ±10.50 to ±30.00 (0.50 steps)	-	-	 (5%, 10%, 20%, 40%)
7.9 to 9.3 (0.2 steps)	12.5, 13.0, 13.5	±10.00 (0.25 steps) ±10.50 to ±30.00 (0.50 steps)	-	-	 (5%, 10%, 20%, 40%)
8.0 to 9.6 (0.4 steps)	14.0	±10.00 (0.25 steps) ±10.50 to ±20.00 (0.50 steps)	-	-	 (10%, 20%)
8.4, 8.8	14.5	±10.00 (0.25 steps) ±10.50 to ±20.00 (0.50 steps)	-0.50 to -6.00 (0.25 steps)	1° to 180° (1° steps)	 (10%, 20%)
8.1 to 9.6 (0.3 steps)	14.0	±25.00 (0.25 steps)	-0.50 to -6.00 (0.25 steps)	1° to 180° (1° steps)	 (5%, 10%, 20%, 40%)

Further CLPL Product Information

CLPL supplies all system and full prescription contact lenses, soft and RGP, in a large range of designs and materials. We also supply fully bespoke contact lenses to the specifications and designs of your choice.

CLPL offers numerous support services to assist customers' businesses. These include EDI links, electronic ordering, private-labelling, bar-coding and direct distribution.

With fully networked CNC manufacturing equipment, CLPL can supply you with any design of your choice in a vast range of materials.

Please contact us for further details.

Our product range includes:

Soft Lenses

Simplon™ Sphere	SAM® MF
Cert™ Toric	Therapeutic
Paediatric	Aphakic

RGP and Other Lenses

Topflex™ Sphere	Eureka™ Asphere
Back Toric	Front Toric
Compensated Toric	Bi-Toric
SAM® MF	Keratoconus
Xtralens™	Apex
Post Graft	Your Design

Ordering Information

To order simply let us know your chosen lens design and prescription (including chosen addition where applicable) or alternatively supply us with K readings, Spectacle Rx, BVD and HVID and we can design the lenses for you. All of CLPL's system lenses carry a full 90-day fitting guarantee.

CLPL accepts orders in a variety of ways including,

Tel: +44 (0) 1223 833936
Fax: +44 (0) 1223 837415
Email: order@clpl.co.uk
Web: www.clpl.co.uk
Other: Contact Lens Calculator and Ordering Program
CLPL EDI systems
Courier and post (see address below)

Visit www.clpl.co.uk or contact us for further information.

We really value your business. We are here to help you!



Contact Lens
Precision Laboratories Ltd

Contact Lens Precision Laboratories Ltd
Mill Lane, Sawston, Cambridge
CB2 4HZ, United Kingdom
Tel: +44 (0)1223 833936 Fax: +44 (0)1223 837415
E-mail: order@clpl.co.uk Web: www.clpl.co.uk