



SILICONE
HYDROGEL
High Myopic™

Now supplied with blue handling tint and UV protection

A 3-monthly contact lens, giving optimum performance and comfort, in a high water content Silicone Hydrogel material for the High Myope.

Healthy & Comfortable

Lathe-cut Silicone Hydrogel material
DK of 60 and very high water retention properties
Provides a more healthy and comfortable lens

Excellent Vision

Designed specifically for the optics of the high myopic eye
Enhanced visual acuity and comfort

Advanced Design

Lens thickness controlled through innovative design
Optimising front optic diameter
Prescription related optical system

Extensive Power Range

Available in powers up to -45.00DS

Now available as a cost effective 2-pack option

Material	26% Filcon V 3, 74% Water
Base Curve	7.00mm to 9.80mm (0.10 steps)
Diameter	12.50mm to 16.00mm* (0.10 steps)
Lens Design	Reduced FOD providing thinner lens
Power Range (D)	-0.25DS to -45.00DS (0.25 steps)
Modulus	0.27 MPa (Typical of mid water content materials)
Handling Tint	Clear
DK	60 x 10 ⁻¹¹ (cm ² /sec) [ml O ₂ /(ml x mm Hg)]
Wear Modality	3-monthly replacement: Daily wear
Pack Size	Single lens, 2-Pack

* Extended range available on request

