



Contact Lens for Keratoconus & all Irregular Corneas



### Life-Changing

Comfortable wear and excellent vision all day, 7 days a week  
Search KeraSoft at [www.youtube.com](http://www.youtube.com) to see Patient Testimonials

### Unique Design

A complete system enabling healthy management of irregular corneas with a completely customisable soft contact lens

### Healthy & Comfortable

Easy to handle 74% water content Silicone Hydrogel 3-month replacement material for increased oxygen permeability and comfort  
Also available in 77% high water content Hydrogel for annual replacement  
Award-winning lens design and optics to provide enhanced visual acuity

### Simple Fitting System

8 lens fitting set supplied in 77% High Water Hydrogel  
New fitting collateral to assist in achieving optimum fit, comfort and vision

### Wide Range of Parameters

Sphere and Toric  
Base Curve: 7.40mm to 9.40mm  
Diameter: 14.00mm to 15.50mm  
Power: +/-30.00DS  
Cyl up to -15.00DC / Axis 1° to 180°  
Periphery options: Standard, Steep 1,2\*,3\*,4\* (\*Reverse Geometry), Flat 1,2,3,4  
Sector Management: Periphery can be individually customised to allow steeper and flatter sectors with a continuous blend between the two sectors

### Material Comparison Table

	Lathe-cut Silicone Hydrogel (Filcon V 3)	UltraVision 77% (Filcon II 2)
<b>Water Content</b>	74%	77%
<b>Dk</b>	60 x 10 <sup>-11</sup> (cm <sup>2</sup> /sec)[mL O <sub>2</sub> /(mL x mmHg)]	53 x 10 <sup>-11</sup> (cm <sup>2</sup> /sec)[mL O <sub>2</sub> /(mL x mmHg)]
<b>Handling Tint</b>	Clear	Clear
<b>Wear Modality</b>	3 month, daily wear	12 month, daily wear

