













ULTRAVISION INNOVATION IN PRACTICE

PREFERRED CHOICE SPECIALIST CONTACT LENSES

Condition	Brand	Material	Modality *	Parameter Range
Keratoconus (Early-Moderate-Advanced) Post Graft Post Refractive Surgery Post Cross Linking	 	Hydrogel 77% Water Silicone Hydrogel 74% Water	12 Month 3 Month	BCOR 7.40 to 9.40mm / TD 14.00 to 15.50mm / +/- 30.00D / Cyl to -15.00D Peripheral Fit FLT1/FLT2/FLT3/FLT4 / STP1/STP2*/STP3*/STP4* (*Reverse Geometry) Sector Management Control (SMC) allows different peripheral fits at any sector around the lens
Astigmatic/Myopia	 High Myopic SiH™ DuraWave® Eureka™  	Silicone Hydrogel 74% Water Silicone Hydrogel 74% Water Hioxifilcon B II 49% Water Rigid - Various	3 Month 3 Month 12 Month 12 Month	BCOR 7.00 to 9.00mm / TD 12.50 to 16.50mm / BVP +/- 30.00D / Cyl to -11.00D BCOR 7.00 to 9.80mm / TD to Rx / BVP up to -45.00D BCOR 7.80 to 9.80mm / TD 13.00 to 15.50mm / BVP +/- 30.00D / Cyl to -11.00D BCOR 7.00 to 9.00mm / TD 8.50 to 10.00mm / BVP +/- 30.00D / Cyl to -6.00D
Paediatric/Aphakic	Paediatric SiH™ Simplon™ Aphakic SiH™ Eureka™ Asphere  	Silicone Hydrogel 74% Water MMA/VP 72% Water Silicone Hydrogel 74% Water Rigid - Various	3 Month 12 Month 3 Month 12 Month	BCOR 6.50 to 8.20mm / TD to Rx / BVP +/- 45.00D BCOR to Rx / TD to Rx / BVP +/- 45.00D BCOR 7.00 to 9.80mm / TD to Rx / BVP up to +45.00D BCOR 7.00 to 9.00mm / TD 8.50 to 10.00mm / BVP +/- 30.00D / Cyl to -6.00D
Epithelial Defects (Promotes Healing)	Bandage SiH™ UltraVision Bandage 77% UltraVision Bandage 67%  	Silicone Hydrogel 74% Water Filcon II 2 77% Water Filcon II 2 67% Water Silicone Hydrogel 58% Water	3 Month 12 Month 12 Month 1 Month	BCOR 7.00 to 9.80mm / TD 12.50 to 19.00mm / BVP Plano BCOR 7.00 to 9.80mm / TD 12.50 to 19.00mm / BVP Plano BCOR 7.00 to 9.80mm / TD 12.50 to 19.00mm / BVP Plano BCOR 8.40mm / TD 14.10mm / BVP +8.00 to -10.00D
Presbyopia	 Multifocal  Multifocal SAM® Multifocal - Soft SAM® Multifocal - GP 	Silicone Hydrogel 74% Water Hioxifilcon B II 49% Water GM3 58% Water Rigid - Various	3 Month 12 Month 12 Month 12 Month	BCOR 7.00 to 9.00mm / TD 12.50 to 16.50mm / BVP +/- 30.00D / Cyl to -11.00D / Add up to +3.00 BCOR 7.80 to 9.80mm / TD 13.00 to 15.50mm / BVP +/- 30.00D / Cyl to -11.00D / Add up to +3.00 BCOR 7.00 to 9.800mm / TD 12.50 to 16.00mm / BVP +/- 30.00D / Cyl to -11.00D / Add to +3.50 BCOR 7.00 to 9.00mm / TD 8.70 to 10.20mm / BVP +/- 30.00D / Cyl to -6.00D / Add to +4.50

* All Lenses are Daily Wear only

Full Rx range available, tailor-made to your design

To order **FREEPHONE 0800 585115**
or visit **www.ultravision.co.uk**

INNOVATION IN PRACTICE

UltraVision Contact and Ordering Information

Orders

Orders can be placed by email, fax, phone, mail or online at www.ultravisionorders.co.uk. Please contact Customer Services for login details.

Customer Services

Our Customer Services team and Clinical Advisors are available: Monday - Friday 8.30am to 6.00pm.

Encore™

Practitioners who register patients on our Encore™ frequent replacement scheme benefit from preferential pricing. A £10 fee will be applied to Encore™ cancellations where less than two weeks notification is given prior to date of recall. There is no cancellation fee in relation to our disposable lenses.

Tinting

Other manufacturers' lenses are tinted at customers' own risk. We do not credit or exchange tinted, coloured or prosthetic lenses, if no fitting warranty has been purchased, unless there is a manufacturing fault.

Delivery

Disposable lenses are despatched within 48 hours of receipt of order. Every effort is made to ensure that all other lenses are despatched within 5 working days from receipt of the order. In some cases, such as tinted lenses, extra time may be required. A minimum charge of £2.50 is made to cover postage and packing.

Return of Product

Unopened stock lenses may be returned for exchange or credit, providing they are in a good saleable condition with a minimum of 2 years shelf life remaining.

Exchange Guarantee

Unless already included, an optional fitting warranty may be purchased upon request excluding disposables and multi-packs. The fees for the warranty are quoted in the price list. The fitting warranty allows for up to two exchanges within 90 days from the original order date. Should the lenses prove unsatisfactory a 75% credit will be given on the total lens cost (excluding the warranty charge). We are unable to credit lenses that are returned due to the patient not collecting their lenses.

Payment

Our standard payment terms are 30 days from the date of the invoice. Failure to observe these terms will lead to the account being placed on hold. The Company reserves the right to cease trading with any customer if there is an unsatisfactory observance of our payment terms or any

Terms and Conditions

These terms and conditions shall be construed in accordance with English law and any dispute shall be subject to the jurisdiction of the English Courts. All orders are deemed to have been placed by the customer and accepted by UltraVision upon and subject to these terms and conditions. All goods supplied remain the property of UltraVision until payment is received in full. Products are supplied in the UK on the basis that they will be sold through recognised optical outlets. UltraVision reserves the right to alter any prices and terms without prior notice. UltraVision also reserves the right to refuse to accept or to cancel any order at any time.

Our Product Range also includes:

Soft Disposable - clear1-day® - Specialty 55AB - Specialty 55 - Specialty Colours

Soft Rx - Simplon™ - Omega 38 - CD38 - Presto 58 - Prima 67 - Igel™ Toric - Igel™ Sphere - MiNi™ - H67™ - RxT™ - Simplon™ Cert - Simplon™ Therapeutic - Simplon™ Aphakic

Soft Multifocal - SAM® MF

Tints - Omniflex Softints™ - Zero6 Softints™ - MiNi Softints™ - H67 Softints™ - Prosthetics

Please refer to the full product sheet for detailed information on any products included on this chart

ULTRAVISION 
INNOVATION IN PRACTICE



UltraVision International Ltd

Commerce Way, Leighton Buzzard, Bedfordshire, LU7 4RW, United Kingdom

UK FREEPHONE 0800 585115

Tel: +44 (0)1525 381112 Fax: +44 (0)1525 370091

Email: orders@ultravision.co.uk

www.ultravision.co.uk

www.siliconehydrogel.co.uk

www.kerasoft.co.uk

www.ultravisionorders.co.uk